



WishBone Carbide Drill Instructions for Use

General - Intended Use

The WishBone Carbide Drill may be used to remove a broken or damaged screw, or may be used in conjunction with the WishBone Medical Screw Extraction Kit.

- **Contents:**
The package contains one surgical instrument for use with the removal of broken or damaged screws.
- **Material:**
The instrument is made from a carbide alloy.
- **MRI Compatibility:**
This instrument has not been evaluated for safety and compatibility in the MR environment and has not been tested for heating or migration in the MR environment.
- **Safety Precautions:**
 - Prior to use, thoroughly read this instruction for use. Keep the instructions for use accessible to all staff.
 - The use of a surgical instrument for tasks other than those for which they are intended may result in damaged or broken instruments or patient injury.
- **Adverse Reactions:**
Possible reactions may include but are not limited to:
 - Clinical failure due to inappropriate usage.
 - Necrosis due to thermal load (power driven tools).

Sterile:   

Sterilized with ethylene oxide gas. Caution: For one procedure only. Do not re-sterilize. Do not use if package is open or damaged. This is a single use device. Re-use of this device can result in the transfer of materials not limited to bone, tissue, blood or infectious disease. This device is provided sterile and re-sterilization of the device has not been validated.

How Supplied – Storage:

This Instrument is packed in protective packaging that is labeled to its contents. All instruments are supplied sterile.

- Always store the instrument in the original protective packaging.
- Store the instruments in a dry and dust-free place (standard hospital environment).

Inspection:

Before use, inspect the instrument box carefully. Do not use when sterile barrier is visually damaged.

Warning:

Please note that WishBone Medical supplies instruments as sterile, single use devices.

Once a single use device (SUD) comes into contact with human blood or tissue reprocessing is required. If the end user wishes to re-use the device(s), re-processing can only be performed by an authorized third-party re-processor who has received a 510(k) clearance for such.

Distributed By:
WishBone Medical, Inc.
2150 North Pointe Drive
Warsaw, IN 46582
(574) 306-4006 (P)
www.WishBoneMedical.com

Manufactured By:
LH Industries
4420 Clubview Drive
Fort Wayne, IN 46804
260-32-5563 (P)
www.lhindustries.com